GenoQuick® MTB

VER 1.0

Deutsch: S. 2-12

English: p. 13-23

Français : p. 24-34

Italiano: p. 35-45

Español: p. 46-56

Português: p. 57-67

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Art.-Nr.: 51610

GenoQuick® MTB

Rapid Molecular Genetic Assay for the Direct Detection of the *Mycobacterium tuberculosis* Complex from Patient Specimens

Methodology

The **GenoQuick® MTB** test permits the fast molecular genetic identification of the *M. tuberculosis* complex directly from decontaminated patient specimens like sputum (induction or expectoration), bronchial material (bronchoalveolar lavages and bronchial aspirates), tracheal secretion, gastric juice, and urine. The amplification product is qualitatively detected on a dipstick. First, the single stranded amplicons hybridize with specific probes included in the Primer Nucleotide Mix. These complexes then selectively bind to the test band on the dipstick and are visualized by a gold labeling. The **GenoQuick® MTB** assay is a direct test for a screening of patients with a possible tuberculosis infection.

The whole procedure is divided into three steps: (i) DNA extraction from decontaminated patient specimens (the necessary **GenoLyse®** kit is not provided; see chapter Ordering Information), (ii) multiplex amplification with differently labeled primers (the necessary thermostable DNA polymerase is not provided) with subsequent hybridization, and (iii) detection on the dipstick.

Storage and Precautions

Store Primer Nucleotide Mix (PNM) and Positive Control DNA (C+ GQ MTB) at 2-8°C upon arrival isolated from any potential source of contaminating DNA. If longer storage (more than 4 weeks) is required, store at -20°C. In order to avoid repeated freezing and thawing, aliquot PNM. Store all other kit components at 2-8°C. Close tube with dipsticks properly after each use; the inner tube wall is coated with desiccant. Do not use the reagents beyond their expiry date. Do not mix reagents from kits with different lot numbers.

Patient specimens must always be considered as potentially infectious. Samples from risk patients must always be labeled and handled under suitable safety conditions. Observe all federal, state, and local safety and environmental regulations. Always wear suitable protective clothing and gloves.

Observe the usual precautions for amplification set-up. It is essential that all reagents and materials used for DNA extraction and amplification set-up are free from DNases.

Quality Control

In order to validate the correct performance of the test and the proper functioning of kit components, each dipstick includes 2 control zones:

- a Conjugate Control zone to check the binding of the conjugate on the dipstick
- an Amplification Control zone to check for a successful amplification reaction

DNA Extraction

Decontaminated patient specimens like sputum (induction or expectoration), bronchial material (bronchoalveolar lavages and bronchial aspirates), tracheal secretion, gastric juice (from adults only), and urine are used as starting material for DNA extraction. The assay must not be used as a standard method for detection of *M. tuberculosis* complex from cultured material. Samples must be processed using the NALC/NaOH method according to the CDC publication "Public Health Mycobacteriology: A Guide for the Level III Laboratory". The working area must be free from amplified DNA.

For DNA extraction the **GenoLyse**[®] kit (DNA extraction via alkaline lysis; see chapter Ordering Information) is used and processed as described in the respective manual. The performance evaluation of the **GenoQuick**[®] **MTB** kit was carried out using the method described above. The user is in charge of testing the applicability of any alternative extraction method with this assay.

Amplification

Prepare the amplification mix (45 µl) in a DNA-free room. The DNA sample should be added in a separate area. **It is crucial to use a hot start DNA polymerase.**

Per tube mix:

- 35 µl PNM – provided
- 5 µl 10x polymerase incubation buffer – not provided
- x µl MgCl₂ solution¹¹ – not provided
- 1-2 unit(s) thermostable DNA polymerase (refer to manual) – not provided
- y µl water to obtain a volume of 45 µl (not considering volume of enzyme) – not provided
- Add 5 µl DNA solution leading to a final volume of 50 µl (not considering volume of enzyme).

¹¹ Depending on the enzyme/buffer system used, the optimal MgCl₂ concentration may vary between 1.5 and 4 mM. Please note that some incubation buffers already contain MgCl₂.

The performance evaluation of the **GenoQuick® MTB** assay was carried out using the HotStarTaq DNA Polymerase from Qiagen. When using this enzyme, the following amounts are required per sample:

- 35 µl PNM – provided
- 5 µl 10x PCR Buffer for HotStarTaq (contains 15 mM MgCl₂) – not provided
- 5 µl 25 mM MgCl₂ solution – not provided
- 0.2 µl (1 U) HotStarTaq – not provided
- 5 µl DNA solution (add in a separated area)

The final MgCl₂ concentration in this amplification mix is 4 mM.

Determine the number of samples to be amplified (number of samples to be analyzed plus control samples). A contamination control sample contains 5 µl of water instead of DNA solution. Prepare a master mix containing all reagents except for DNA solution and mix well (do not vortex). Aliquot 45 µl in each of the prepared PCR tubes.

Amplification profile²⁾:

15 min 95°C 1 cycle

30 sec 95°C	}	40 cycles
40 sec 55°C		
30 sec 72°C		

2 min 95°C 1 cycle


5 min 20°C 1 cycle

(∞ 8°C)

²⁾ Applies to the Taq polymerase used for validation. When using other hot start DNA polymerases, the time interval of the first step might have to be reduced (please refer to manual of the enzyme).

Amplification products can be stored at +8 to -20°C. Please note, however, that the hybridization product has only a limited stability. Therefore, repeat the hybridization steps of the cycler program (2 min 95°C, 5 min 20°C) if the time interval between amplification and detection exceeds 30 min at room temperature (15-25°C) and 1 d at +8 to -20°C, respectively.

Detection

Note: The dipstick is not displayed in original size.

Preparation

Directly before starting the detection, remove one dipstick from the dipstick tube for each sample and mark it with a pencil in the white label area. The dipstick is almost completely covered with a transparent foil and can therefore be touched at the label area. For each sample, pipette 100 μ l Running Buffer (RB) in a tube of the test rack provided.

- 1. Add 10 μ l hybridized amplicon to a tube filled with 100 μ l RB and pipette up and down to mix well.**

The color of the RB will change to pale yellow after adding the amplicon.

- 2. Place dipstick with the gold matrix facing downwards (arrows pointing upwards) in the respective tube and incubate for 10 minutes at room temperature.**

Positive results may be detected earlier; still, the indicated time should be observed. An extended incubation time, on the other hand, may lead to a change in signal intensity.

3. Remove dipsticks from tubes to read and note down results on the evaluation sheet provided (see chapter Evaluation and Interpretation of Results) and discard dipsticks.

Preserving dipsticks bears a contamination risk; additionally, the drying process may cause changes in signal intensities. Hence, after documenting the results, the developed dipsticks should be discarded or archived digitally.

Dispose of dipsticks and tubes with remaining buffer in a container with freshly prepared 1.5% sodium hypochlorite solution and discard gloves. Treat working area with a freshly prepared 1.5% sodium hypochlorite solution and rinse with water.

Evaluation and Interpretation of Results

An evaluation sheet is provided with the kit.

Conjugate Control (CC)

A line must develop in this zone, documenting the efficiency of conjugate binding and proper flow of the dipstick.

Amplification Control (AC)

When the test is performed correctly, a control amplicon will bind to the Amplification Control zone on the dipstick. If this band is developed, the carry-over of amplification inhibitors and mistakes during setup and performance of the amplification reaction can be excluded.


In case of a positive test result, the signal of the Amplification Control zone can be weak or even vanish totally. This might be due to competition of the single reactions during amplification or the presence of a limited amount of amplification inhibitors. In this case, however, the amplification reaction was performed correctly and the test does not have to be repeated.

In case of a negative test result, the intensity of the AC should be about as strong as the AC of the contamination control sample. A weak or missing AC band in case of a negative test result indicates mistakes during amplification set-up, carry-over of amplification inhibitors, or inhibition due to overload of the test system (see chapter Troubleshooting). In this case, the test result is not valid and the respective sample has to be repeated.

***Mycobacterium tuberculosis* complex (MTB)**


This reaction zone documents the presence of a member of the *M. tuberculosis* complex. The intensity of this band varies depending on the number of cells in the patient sample. In case of a very weak signal, the respective sample should be repeated.

The figure below shows possible valid banding patterns:

Note: The band intensities displayed in the figures above do not provide reference signal strengths for evaluation of test results, but intend to give an impression of the possible variety of band intensities.

The figure below shows possible banding patterns representing invalid results:

not valid

Note: The band intensities displayed in the figures above do not provide reference signal strengths for evaluation of test results, but intend to give an impression of the possible variety of band intensities.

Limitations

The test was validated for usage of decontaminated direct patient material. Clinical specimens tested are sputum (induction or expectoration), bronchial material (bronchoalveolar lavages and bronchial aspirates), tracheal secretion, gastric juice (from adults only), and urine. The user is in charge of validating the applicability of this assay for other sample materials.

As any DNA detection method the test system on hand detects DNA from viable and non-viable bacteria. Therefore, the **GenoQuick® MTB** test may not be used for monitoring the progression or success of treatment of patients with antimicrobial therapy.

Application of exceeded DNA concentrations may lead to inhibition of PCR or detection; therefore, the test must not be used as standard method for examining *Mycobacterium* cultures.

The **GenoQuick® MTB** test generates qualitative results. The intensity of the MTB band on a dipstick does not give information about the number of cells in a sample. A negative test result does not exclude the possibility that the concentration of the target organism lies beyond the detection limit of the assay. As with other diagnostic assays, the results of this direct test may only be interpreted in combination with additional laboratory and clinical data available to the responsible physician.

Prior to amplification, bacterial DNA has to be extracted using the **GenoLyse®** kit (see chapters DNA Extraction and Ordering Information). It must be ensured that the template DNA is efficiently amplified during the amplification reaction.

The test only works within the limits of the genomic regions the primers and probes were chosen from. Potential sequence analysis remains to further investigations. As with any detection system based on hybridization the test system on hand bears the possibility that sequence variations in the genomic regions the primers and probes were chosen from but the detection of which the test was not designed for may lead to false results. Due to the high variability of bacterial genomes it is possible that very rare species variations might not be detected. The test reflects the state of knowledge of Hain Lifescience.

Use of this assay is limited to qualified personnel well trained in the test procedure and familiar with molecular biological methods.

Performance evaluation of this assay was carried out using the HotStarTaq DNA Polymerase from Qiagen (Hilden, Germany). For DNA extraction the **GenoLyse®** kit was used. Since the performance characteristics of this assay have not been validated for all polymerases and DNA extraction kits commercially available, the user is in charge of validating the applicability of polymerases or extraction methods other than those mentioned above.

Performance data of the assay can be requested through:
www.hain-lifescience.com

Troubleshooting

Overall weak or no signals (including Conjugate Control zone)

- Too little or too much Running Buffer added to tube.
- Too much amplicon added to Running Buffer.
- Dipstick not correctly dipped into Running Buffer.

Weak or no Amplification Control with no MTB signal present

- Amplicon not denatured and hybridized. Repeat hybridization (see chapter Amplification).
- Time between denaturation and detection longer than 30 minutes or samples were warmed above room temperature. Repeat hybridization (see chapter Amplification).
- No or too little amount of amplicon added to Running Buffer.
- Quality of extracted DNA does not allow an efficient amplification. Repeat DNA extraction and amplification reaction; if necessary, try a different DNA extraction method (see chapter DNA Extraction).
- DNA concentration used for amplification too high e.g. due to the use of DNA from cultured material that has not been sufficiently diluted

Unexpected result

- Contamination of extracted DNA and/or amplification agents with extracted and/or amplified DNA. In case amplification agents are contaminated, a contamination control sample also shows the respective banding pattern.
- Contamination of neighboring tubes by spillage during addition of amplicon.
- No or too little DNA added to amplification reaction, e.g. due to no or too little cells in patient specimen.
- The bacterial species present in the sample can not be detected with this test.

Material Required but not Provided

- Adjustable pipettes for 10, 20, 200, and 1000 μ l
- Disposable gloves
- Disposable sterile pipette tips with filter
- **GenoLyse**[®] kit (see chapter Ordering Information)
- PCR tubes, DNase and RNase free
- Sodium hypochlorite solution (household bleach without dyes and scents)
- Thermal cycler (heating rate: 3°C/sec, cooling rate: 2°C/sec, precision: +/-0.2°C)
- Thermostable DNA polymerase with buffer (recommendation: hot start enzyme, extension rate: 2-4 kb/min at 72°C, half-life: 10 min at 97°C, 60 min at 94°C, amplification efficiency: >10⁵ fold)
- Timer
- Water (molecular biology grade)

Kit Contents

	12	96
Dipsticks coated with specific probes (DS GQ MTB)	12	96
Primer/Nucleotide Mix (PNM GQ MTB) contains specific oligonucleotides, nucleotides, dye	0.5	4 ml
Running Buffer (RB) ready to use contains buffer, <1% NaCl, <1% anionic tenside	2 ml	16 ml
Test rack with 96 tubes, Manual	1 of each	1 of each
Evaluation sheet	1	4

Ordering Information

GenoLyse[®] DNA extraction kit for 96 samples

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Hain Lifescience GmbH