GenoType MTBC

VER 1.X

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IFU-301-09

CE

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Diagnostic performance

DNA Extraction with quick method

Diagnostic performance characteristics were determined in a study with 82 clinical isolates grown on solid medium [6] and 77 clinical isolates grown in liquid medium [7]. DNA extraction was performed using the quick method described in chapter DNA Extraction. The **GenoType MTBC** was compared to various biochemical and molecular differentiation methods (e.g. determination of colony morphology, nitrate reduction on modified Dubos broth, niacin accumulation test, growth in the presence of thiophene-2-carboxylic acid hydrazide (TCH), spoligotyping, and RD1-PCR).

In the first study, 82 clinical isolates grown on **solid medium** were analyzed [6]. For all samples, the **GenoType MTBC** results were in agreement with the results of the biochemical or molecular differentiation (see table 1).

Table 1: Results of the *M. tuberculosis* complex differentiation [bacteria grown on solid medium; DNA extraction with quick method]

result GenoType MTBC	number of isolates	result conventional methods	
<i>M. tuberculosis</i>	31	31 <i>M. tuberculosis</i>	diagnostic sensitivity: 100% diagnostic specificity: 100% positive predictive value: 100% negative predictive value: 100%
<i>M. africanum</i>	4	4 <i>M. africanum</i>	
<i>M. bovis</i> subsp. <i>bovis</i>	17	17 <i>M. bovis</i> subsp. <i>bovis</i>	
<i>M. bovis</i> subsp. <i>caprae</i>	17	17 <i>M. bovis</i> subsp. <i>caprae</i>	
<i>M. bovis</i> BCG	4	4 <i>M. bovis</i> BCG	
<i>M. microti</i>	9	9 <i>M. microti</i>	
total	82	82	

In the second study, 77 clinical isolates grown in **liquid medium** were analyzed [7]. For all samples, the **GenoType MTBC** results were in agreement with the results of the biochemical or molecular differentiation (see table 2).

Table 2: Results of the *M. tuberculosis* complex differentiation [bacteria grown in liquid medium; DNA extraction with quick method]

result GenoType MTBC	number of isolates	result conventional methods	
<i>M. tuberculosis</i>	71	71 <i>M. tuberculosis</i>	diagnostic sensitivity: 100% diagnostic specificity: 100% positive predictive value: 100% negative predictive value: 100%
<i>M. africanum</i>	1	1 <i>M. africanum</i>	
<i>M. bovis</i> subsp. <i>bovis</i>	5	5 <i>M. bovis</i> subsp. <i>bovis</i>	
total	77	77	

The Primer Nucleotide Mix (PNM) has been replaced by new kit constituents, namely **Amplification Mixes A** and **B** (AM-A and AM-B). In order to check if the change of kit constituents impacts test results, 72 liquid and solid cultures were tested with both kit variants. DNA was extracted using the quick method described in chapter DNA Extraction. In addition, 2 negative controls (water instead of DNA added to the aliquoted master mix) were tested. The results were identical for all samples (see table 3).

Table 3: Results of the *M. tuberculosis* complex differentiation [Primer Nucleotide Mix compared to Amplification Mixes; DNA extraction with quick method]

result GenoType MTBC with Primer Nucleotide Mix	number of isolates	result GenoType MTBC with Amplification Mixes
<i>M. tuberculosis</i>	60	60 <i>M. tuberculosis</i>
<i>M. africanum</i>	4	4 <i>M. africanum</i>
<i>M. bovis</i> subsp. <i>bovis</i>	4	4 <i>M. bovis</i> subsp. <i>bovis</i>
<i>M. bovis</i> subsp. <i>caprae</i>	1	1 <i>M. bovis</i> subsp. <i>caprae</i>
<i>M. bovis</i> BCG	2	2 <i>M. bovis</i> BCG
high GC gram-positive bacterium	3	3 high GC gram-positive bacterium
negative	2	2 negative
total	74	74

DNA Extraction with **GenoLyse**®

In a study comprising 52 *Mycobacterium*-positive cultures (growth on Loewenstein-Jensen medium or in BACTEC MGIT 960 medium), DNA was extracted with the **GenoLyse**® kit and then tested with the **GenoType MTBC**. For comparison, DNA was isolated with the quick method in parallel and then tested with the **GenoType MTBC**.

With both extraction methods, identical results were obtained (see table 4).

Table 4: Results of the *M. tuberculosis* complex differentiation [DNA extraction with quick method compared to DNA extraction with **GenoLyse**®]

result after DNA extraction with quick method	number of isolates	result after DNA extraction with GenoLyse ®
<i>M. tuberculosis</i>	47	47 <i>M. tuberculosis</i>
<i>M. bovis</i> BCG	2	2 <i>M. bovis</i> BCG
high GC gram-positive bacterium	3	3 high GC gram-positive bacterium
total	52	52

The Primer Nucleotide Mix (PNM) has been replaced by new kit constituents, namely **Amplification Mixes A** and **B** (AM-A and AM-B). In order to check if the change of kit constituents impacts test results, 92 solid cultures were tested with both kit variants. DNA was extracted using the **GenoLyse**® kit. In addition, 2 negative controls (water instead of DNA added to the aliquoted master mix) were tested. The results were identical for all samples (see table 5).

Table 5: Results of the *M. tuberculosis* complex differentiation (Primer Nucleotide Mix compared to Amplification Mixes; DNA extraction with **GenoLyse**[®])

result GenoType MTBC with Primer Nucleotide Mix	number of isolates	result GenoType MTBC with Amplification Mixes
<i>M. tuberculosis</i>	48	48 <i>M. tuberculosis</i>
<i>M. bovis</i> subsp. <i>bovis</i>	2	2 <i>M. bovis</i> subsp. <i>bovis</i>
<i>M. bovis</i> BCG	2	2 <i>M. bovis</i> BCG
high GC gram-positive bacterium	40	40 high GC gram-positive bacterium
negative	2	2 negative
total	94	94

Further performance studies of the **GenoType MTBC** were published [8,9,10].

Analytical performance

Analytical specificity

The specificity of the **GenoType MTBC** test is ensured by the accurate design of specific primers and probes which considers, among others, homology comparisons of the sequences published in gene databases, and by stringent reaction conditions.

The analytical specificity of the **GenoType MTBC** (with Primer Nucleotide Mix) was determined with 98 strains of 75 different species and subspecies. The study included the *Mycobacterium tuberculosis* complex strains *M. africanum*, *M. bovis* BCG (5x), *M. bovis* subsp. *bovis* (7x), *M. bovis* subsp. *caprae*, *M. microti*, and *M. tuberculosis* (2x) as well as the following strains not detectable with the test system: *Actinomyces naeslundii*, *Aggregatibacter actinomycetemcomitans*, *Bacillus cereus*, *Bacteroides forsythus*, *Corynebacterium ammoniagenes*, *C. bovis*, *C. durum*, *Escherichia coli*, *Enterococcus faecalis*, *Gordona rubropertinctus*, *Klebsiella oxytoca*, *M. abscessus*, *M. alvei* (2x), *M. asiaticum*, *M. avium*, *M. celatum*, *M. chelonae*, *M. duvalii*, *M. feldmannii*, *M. fortuitum* (2 sequevars), *M. gastri*, *M. genavense*, *M. goodii*, *M. gordonae*, *M. haemophilum*, *M. heckeshornense*, *M. immunogenum* (2x), *M. interjectum* (2x), *M. intermedium* (2x), *M. intracellulare* (2x), *M. kansasii*, *M. lentiflavum*, *M. mageritense*, *M. malmoense*, *M. marinum*, *M. mucogenicum*, *M. nonchromogenicum*, *M. palustre* (2x), *M. peregrinum*, *M. phlei*, *M. ratisbonense*, *M. scrofulaceum*, *M. shimoidei*, *M. simiae*, *M. smegmatis*, *M. szulgai*, *M. terrae*, *M. triplex* (2x), *M. ulcerans*, *M. vaccae*, *M. xenopi*, MRSA, *Nocardia abscessus*, *N. africana*, *N. amarae*, *N. asteroides*, *N. farcinica*, *Porphyromonas gingivalis*, *Prevotella intermedia*, *Rhodococcus equi* (3x), *R. erythropolis*, *Saccharomonospora glauca*, *Staphylococcus aureus*, *S. epidermidis*, *S. intermedium*, *Tannerella forsythia*, *Treponema denticola*, *Tsukamurella inchonensis*, *T. paurometabola* (3x).

The 17 *M. tuberculosis* complex strains were correctly identified. The further 81 isolates displayed no MTBC band and no species-specific banding pattern. Hence, the analytical specificity was 100%.

Reproducibility

Intra-assay precision

In order to determine the intra-assay precision of the **GenoType MTBC** (with Primer Nucleotide Mix), one *M. microti* strain, one *M. bovis* BCG strain, one *M. gordonae* strain, and one negative control were set up in five parallels and tested under identical conditions. DNA extraction was performed using the quick method described in chapter DNA Extraction. All strains showed the expected signals in all five parallels and the negative controls were negative. Hence, the intra-assay precision was 100%.

In order to determine the intra-assay precision of the **GenoType MTBC** (with Amplification Mixes), two *M. bovis* BCG strains, one *M. gordonae* strain, and one negative control were set up in four parallels and tested under identical conditions. DNA extraction was performed using the **GenoLyse**[®] kit. All strains showed the expected signals in all four parallels and the negative controls were negative. Hence, the intra-assay precision was 100%.

Inter-assay precision

In order to determine the inter-assay precision of the **GenoType MTBC** (with Primer Nucleotide Mix), one *M. microti* strain, one *M. bovis* BCG strain, one *M. gordonae* strain, and one negative control were set up in five parallels and tested at three different points in time. The other experimental conditions (instruments, lot numbers, operator, etc.) were identical. DNA extraction was performed using the quick method described in chapter DNA Extraction. All strains showed the expected signals in all five parallels and the negative controls were negative. Hence, the inter-assay precision was 100%.

In order to determine the inter-assay precision of the **GenoType MTBC** (with Amplification Mixes), two *M. bovis* BCG strains, one *M. gordonae* strain, and one negative control were tested at three different points in time. The other experimental conditions (instruments, lot numbers, operator, etc.) were identical. DNA extraction was performed using the **GenoLyse**[®] kit. All strains showed the expected signals and the negative controls were negative. Hence, the inter-assay precision was 100%.

Interfering substances

There are substances that may inhibit PCR reactions. Such inhibitors may, for example, originate from the culture medium. In order to assess if the medium influences the **GenoType MTBC**, 5 different *M. tuberculosis* samples and 5 different *M. bovis* subsp. *bovis* samples were cultured in 4 different media (solid media: Loewenstein-Jensen, Stonebrink, and Middlebrook-7H10, liquid medium: MGIT (Becton Dickinson)). Then the culture samples and 5 non-inoculated media samples were tested with the **GenoType MTBC** (with Primer Nucleotide Mix). All *M. tuberculosis* samples and all *M. bovis* subsp. *bovis* samples showed the correct results and the non-inoculated controls were negative. Hence, it can be excluded that the tested media import inhibitors into the **GenoType MTBC** test.

Other inhibitors may be imported by reagents used for DNA extraction. When using the recommended quick protocol (see chapter DNA Extraction), however, this is impossible because the method is based on mechanical cell lysis where no reagents are added.

Stability

Stability is determined according to DIN EN ISO 23640.

Shelf life of the **GenoType MTBC** test kit when stored as recommended: see box label.

References

1. World Health Organization. Global tuberculosis control: WHO report 2011. WHO/HTM/TB/2011.16. World Health Organization, Geneva, Switzerland 2011.
2. Biosafety in microbiological and biomedical laboratories, 5th edition. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 2009.
3. Protection of laboratory workers from instrument biohazards and infectious disease transmitted by blood, body fluids, and tissue. Approved guideline. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards), USA, Document M29 (please refer to the latest version).
4. Kent PT, Kubica GP. Public health mycobacteriology: a guide for the level III laboratory. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 1985.
5. Isenberg HD. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C., USA 1992.
6. Richter E, Weizenegger M, Rüschi-Gerdes S, Niemann S. Evaluation of GenoType MTBC assay for differentiation of clinical *Mycobacterium tuberculosis* complex isolates. *J Clin Microbiol* 2003; 41:2672-2675.
7. Richter E, Weizenegger M, Fahr AM, Rüschi-Gerdes S. Usefulness of the GenoType MTBC assay for differentiating species of the *Mycobacterium tuberculosis* complex in cultures obtained from clinical specimens. *J Clin Microbiol* 2004; 42:4303-4306.
8. Gomez MP, Herrera-León L, Jiménez MS, Rodríguez JG. Comparison of GenoType MTBC with RFLP-PCR and multiplex PCR to identify *Mycobacterium tuberculosis* complex species. *Eur J Clin Microbiol Infect Dis* 2007; 26:63-66.
9. Neonakis IK, Gitti Z, Petinaki E, Maraki S, Spandidos DA. Evaluation of the GenoType MTBC assay for differentiating 120 clinical *Mycobacterium tuberculosis* complex isolates. *Eur J Clin Microbiol Infect Dis* 2007; 26:151-152.
10. Somoskovi A, Dormandy J, Rivenburg J, Pedrosa M, McBride M, Salfinger M. Direct comparison of the GenoType MTBC and genomic deletion assays in terms of ability to distinguish between members of the *Mycobacterium tuberculosis* complex in clinical isolates and in clinical specimens. *J Clin Microbiol* 2008; 46:1854-1857.

Important Changes in IFU-301-09

Chapter	Change
	generally revised and restructured: <ul style="list-style-type: none">– former chapter "Methodology" extended and split in new chapters "Intended Use" and "Principles of the Procedure"– former chapter "Storage and Precautions" extended and split in new chapters "Storage and Disposal of Kit Constituents", "Precautions for Handling Kit Constituents", and "Specimen Requirements"– new chapters: "Summary and Explanation", "Ordering Information", "Performance Characteristics", "References", "Important Changes"
DNA Extraction several (especially Amplification)	DNA extraction can now be done using the GenoLyse® kit. In order to obtain a fully reactive master mix, now only two pipetting steps have to be performed combining the two Amplification Mixes A and B (AM-A and AM-B). They come as Kit Component 2 and must be stored at -20°C. AM-A already includes Taq polymerase and buffer; hence, the enzyme has not to be purchased and added separately any more.

HAIN

LIFESCIENCE

301-09-03

CE

IVD

Hain Lifescience GmbH

Hardwiesenstraße 1, 72147 Nehren, Germany
www.hain-lifescience.de, +49 (0) 74 73- 94 51- 0