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
1. APPELBAUM P.C., LEATHERS D.J.
Evaluation of the Rapid NFT System for Identification of Gram-Negative, Nonfermenting Rods.
(1984) J. Clin. Microbiol. 20, 730-734.
2. BERNARDS A.T., VAN DER TOORN J., VAN BOVEN C.P.A., DIJKSHOORN L.
Evaluation of the Ability of a Commercial System to Identify *Acinetobacter* Genomic Species.
(1996) Eur. J. Clin. Microbiol. Infect. Dis. 15, 4, 303-308.
3. BILKEY M.K., BREMNER D.A., CAMERON G.L., GARNER J.G.
Comparison of Five Commercial Methods for the Identification of Non-fermentative and Oxidase Positive Fermentative Gram-Negative Bacilli.
(1988) N.Z.J. Med. Lab. Technol., 8-12.
4. CULLEN K.C., KLOOSTERMAN R.E., SHALIS P.J., PIERSON C.L.
Comparison of Three Commercial Systems for the Identification of Glucose Non-fermenting Gram-Negative Rods.
(1989) ASM Annual Meeting - Poster N° C26.
5. DANCE D.A.B, WUTHIEKANUN V., NAIGOWIT P., WHITE N.J.
Identification of *Pseudomonas pseudomallei* in Clinical Practice : Use of Simple Screening Tests and API 20 NE.
(1989) J. Clin. Pathol. 42, 645-648.
6. FRENEY J., GARONNAT D., BOUVARD V., FLEURETTE J.
Comparaison de Deux Systèmes d'Identification de Bacilles Gram Négatifs Non Fermentants et de Bacilles Fermentants Oxydase Positive.
(1984) Ann. Biol. Clin. 42, 337-341.
7. GEISS H.K., PIOTROWSKI H.D., HINGST V.
Evaluation of API 20 NE in Routine Diagnostics of Nonfermenting Gram-Negative Rod-Shaped Bacteria.
(1985) Zbl. Bakt. Mikrobiol. Hyg. A-Med. 259, 1, 35-42.
8. KISKA D.L., KERR A., JONES M.C., CARACCILO J.A., ESKRIDGE B., JORDAN M., MILLER S., HUGHES D., KING N., GILLIGAN P.H.
Accuracy of Four Commercial Systems for Identification of *Burkholderia cepacia* and Other Gram-Negative Nonfermenting Bacilli recovered from Patients with Cystic Fibrosis.
(1996) J. Clin. Microbiol. 34, 4, 886-891.
9. LAMPE A.S., VAN DER REIJDEN T.J.K.
Evaluation of Commercial Test Systems for the Identification of Nonfermenters.
(1984) Eur. J. Clin. Microbiol. 3, 301-305.
10. MARTIN R., SIAVOSHI F., McDOUGAL D.L.
Comparison of Rapid NFT System and Conventional Methods for Identification of Nonsaccharolytic Gram-Negative Bacteria.
(1986) J. Clin. Microbiol. 24, 1089-1092.
11. OVERMAN T.L., KESSLER J.F., SEABOLT J.P.
Comparison of API 20 E, API Rapid E and API Rapid NFT for Identification of Members of the Family *Vibrionaceae*.
(1985) J. Clin. Microbiol. 22, 778-781.
12. PALMIERI M.J., CARITO S.L., MEYER R.F.
Comparison of Rapid NFT and API 20 E with Conventional Methods for Identification of Gram-Negative Nonfermentative Bacilli from Pharmaceuticals and Cosmetics.
(1988) Applied and Environmental Microbiol. 54, 2838-2841.
13. PELADAN F., MONTEIL H.
Identification of *Pseudomonas*, *Flavobacterium* and *Alcaligenes* with the API 20 NE System.
(1988) Path. Biol., 36, 187-192.
14. SOGAARD P., GAHRN-HANSEN B., HUI-PING Z., FREDERIKSEN W.
An Investigation of Three Commercial Methods for Rapid Identification of Non-Enteric Gram-Negative Rods.
(1986) Acta Path. Microbiol. Immunol. Scand. Sect. B, 94, 357-363.
15. TOWNER K.J., CHOPADE B.A.
Biotyping of *Acinetobacter calcoaceticus* using the API 20 NE System.
(1987) Journal of Hospital Infection. 10, 145-151.
16. VON GRAEVENITZ A., ZOLLINGER-ITEN J.
Evaluation of Pertinent Parameters of a New Identification System for Non-Enteric Gram-Negative Rods.
(1985) Eur. J. Clin. Microbiol. 4, 108-112.
17. Clinical and Laboratory Standards Institute, M50-A Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 n° 23.

TABLE DES SYMBOLES / INDEX OF SYMBOLS / SIMBOLE / CUADRO DE SIMBOLOS /
TABELLA DEI SIMBOLI / QUADRO DOS SÍMBOLOS / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ / SYMBOLER /
SYMBOLFORTEGNELSE / TABELA SYMBOLI

Symbole / Symbol Símbolo / Simbolo Σύμβολο	Signification / Meaning / Bedeutung Significado / Significato / Επεξήγηση Betydelse / Betydning / Znaczenie
 REF	Référence du catalogue Catalogue number (GB) / Catalog number (US) Bestellnummer / Número de catálogo / Numero di catalogo Referência de catálogo / Αριθμός καταλόγου Katalognummer / Katalognummer / Numer katalogowy
	Dispositif médical de diagnostic in vitro In Vitro Diagnostic Medical Device / In Vitro Diagnostikum Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro Dispositivo médico para diagnóstico in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Medicintekniska produkter för in vitro diagnostik Medicinsk udstyr til in vitro-diagnostik Wyrób do diagnostyki In Vitro
	Fabricant / Manufacturer / Hersteller / Fabricante Fabbicante / Κατασκευαστής / Tillverkare / Producent
	Limites de température / Temperature limitation Temperaturbegrenzung / Limite de temperatura Limiti di temperatura / Limites de temperatura Περιορισμοί θερμοκρασίας / Temperaturbegränsning Temperaturbegrænsning Przestrzegać zakresu temperatury
	Utiliser jusque / Use by / Verwendbar bis Fecha de caducidad / Utilizzare entro / Prazo de validade Ημερομηνία λήξης / Använd före / Holdbar til / Użyć przed
	Code du lot / Batch code Chargenbezeichnung / Código de lote Codice del lotto / Código do lote Αριθμός Παρτίδας / Lot nummer / Lotnummer / Kod partii
	Consulter les instructions d'utilisation Consult Instructions for Use Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultare le istruzioni per l'uso Consulte as instruções de utilização Συμβουλευτείτε τις οδηγίες χρήσης Se handhavandebeskrivningen / Se brugsanvisning Sprawdź w instrukcji obsługi
	Contenu suffisant pour "n" tests Contains sufficient for <n> tests Inhalt ausreichend für <n> Prüfungen Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Conteúdo suficiente para "n" ensaios Περιεχόμενο επαρκές για «n» εξετάσεις Räcker till "n" antal tester Indeholder tilstrækkeligt til "n" test Wystarczy na wykonanie <n> testów